

## SER: Seminal Fluid – Acid Phosphatase Spot Test

---

**Principle** Acid phosphatase (AP) is an enzyme that is found in seminal fluid, a component of semen. Its concentration in seminal fluid is up to 400 times greater than that found in other body fluids.

Sodium  $\alpha$ -naphthyl phosphate is cleaved by acid phosphatase in the sample, releasing sodium phosphate and naphthol. Naphthol couples with brentamine (*o*-dianisidine) to create a purple azo dye. The formation of a purple color indicates the presence of acid phosphatase.

---

**Supplies** The following supplies are used in this procedure:

- filter paper
  - cotton swabs
  - glass plates
- 

**Reagents** This procedure uses the following reagents:

- AP Spot Test reagent
  - Dissolve 0.26 grams of SERI Acid Phosphatase Spot Test PMR in 10 mL of deionized water. This solution must be made fresh prior to use.

---

**Quality control** The *AP Spot Test reagent* must be quality control tested before each use with

- a positive control (semen or seminal stain)
- a negative control (water)

---

**Records** Record the lot number and expiration date of the reagent, as well as the results of the quality control tests, in the examination documentation.

---

*Continued on next page*

## SER: Seminal Fluid – Acid Phosphatase Spot Test, Continued

### Overlay procedure

Use the following procedure to perform the test using filter paper overlays.

Step	Action
1	Moisten a suitably sized piece of filter paper with deionized water.
2	Lay the filter paper over the suspected stain to transfer it to the filter paper.  A glass or plastic plate and weight can be applied to ensure proper contact.  Mark seams or other reference points for orientation.
3	After removing the filter paper from the evidence item, apply the <i>AP Spot Test reagent</i> until the filter paper is saturated.
4	Record the time of development of any color that occurs within 60 seconds.  If no color reaction occurs within 60 seconds, record the result as negative.

### Procedure using swabs or filter paper

Use the following procedure to perform the test using swabs or filter paper.

Step	Action
1	Moisten a sterile test swab or filter paper with deionized water.
2	Apply the test swab or filter paper directly onto the item to be tested.
3	Add a drop of <i>AP Spot Test reagent</i> to the test swab or filter paper.
4	Record the time of development of any color that occurs within 60 seconds.  If no color reaction occurs within 60 seconds, record the result as negative.

*Continued on next page*

## SER: Seminal Fluid – Acid Phosphatase Spot Test, Continued

---

**Interpretation**      The development of a purple color within 60 seconds is a positive (+) result for the presence of acid phosphatase.

The absence of a color reaction within 60 seconds is a negative (-) result for the presence of acid phosphatase.

Acid phosphatase is present in seminal fluid, a component of semen, in high concentrations; however, it is also present in other body fluids at lower concentrations and in plants, fungi, and bacteria. Semen stains tend to give a faster and stronger reaction than other sources.

Since acid phosphatase is not unique to seminal fluid, the test is only a presumptive test for the presence of seminal fluid.

A negative reaction does not necessarily mean that other components of semen are not present. If there is a negative reaction for acid phosphatase and no spermatozoa are detected in the sample extract, further testing for p30 or semenogelin must be conducted.

---