

SER: Preservation of Liquid Blood

Principle Liquid blood samples are preserved as dried stains on bloodstain cards.

Equipment and supplies This procedure uses the following laboratory equipment and supplies:

- Whatman® FTA bloodstain cards
 - disposable pipettes
 - manila envelopes.
-

Procedure Use the following procedure to preserve liquid blood samples.

Step	Action
1	Transfer the liquid blood to a bloodstain card. NOTE: The blood card is labeled with the appropriate information including the laboratory case number and item number.
2	When the blood is completely dry, place in a manila envelope, seal envelope, and return to kit.
