

## SER: Species Determination- Crossover Electrophoresis

---

### Principle

Crossover electrophoresis is used to determine the species origin of body fluids and tissues.

In this variation of the precipitin method, antigen and precipitating antiserums are brought together electrophoretically.

- The antigens in the fluid or tissue extracts are serum albumin and  $\alpha$  and  $\beta$  globulins which, under the test conditions, move toward the anode.
  - The antibodies in the antiserums are  $\gamma$  globulins and are moved by electroendosmosis toward the cathode.
  - Reacting antibody and antigen form a precipitate that can be visualized with a protein stain.
- 

### Safety

Electrophoretic procedures have electrical hazards associated with them. Therefore, appropriate caution should be exercised when implementing this procedure.

---

### Equipment and supplies

This procedure uses the following laboratory equipment and supplies:

- electrophoresis tank
  - electrophoresis power supply
  - analytical balance
  - pH meter
  - oven
  - Bunsen burner, heating plate, or microwave
  - glass or plastic support, approximately 7.5 cm x 5 cm
  - pipette
  - pipette tips
  - microcentrifuge tubes
  - gel bond
  - filter paper
  - weight (optional)
  - stereomicroscope
- 

*Continued on next page*

## SER: Species Determination- Crossover Electrophoresis, Continued

---

### Reagents

This procedure uses the following reagents:

- 20% w/v Sodium Hydroxide
    - Dissolve 20 grams of sodium hydroxide in deionized water to bring to a volume of 100 mL.
  - 1 M NaCl
    - Dissolve 58.4 grams of NaCl in deionized water to bring to a volume of 1 liter.
  - Destain Solution
    - Mix the following: methanol/acetic acid/deionized water (50/10/50).
  - Stain Solution
    - Dissolve 0.2 grams of Brilliant Blue in 100 mL of Destain Solution.
  - human and animal antiserums
  - animal normal serums
- 

### Buffer

This procedure uses the following buffer:

For 1 liter, mix:

- 21.8 g glycine
- 4.5 g Trisma base

Dissolve the chemicals in deionized water and bring to final volume to 1 liter. Check the pH and adjust to 8.4 using 20% w/v Sodium Hydroxide (see *Reagents*, above).

Store buffer solution in refrigerator at  $5^{\circ}\text{C} \pm 5^{\circ}\text{C}$ .

If problems are noted while performing this procedure (for example, no precipitin bands are visible between the human standard and the anti-human serum), discard the buffer.

NOTE: after every five uses, discard the buffer, rinse the tank, and refill the tank with buffer.

---

*Continued on next page*

## SER: Species Determination- Crossover Electrophoresis, Continued

---

**Support  
medium**

Using heat, dissolve 0.10 grams of agarose (SERI type E25 or agarose of similar EEO) in 10 mL of buffer.

---

**Standards and  
controls**

The following standards and controls must be used in each test:

- known human blood or serum (human or other species, depending upon the species being tested for)
    - if using known human blood, dilute it with deionized water to red-brown color
  - extractant blank(s) (deionized water or 5% Ammonium Hydroxide).
- 

**Serum and  
antiserum  
quality control**

Each lot of normal serum and antiserum received by the laboratory must be quality control tested prior to use in casework.

- The normal serum is titrated against the antiserum using a serial dilution of normal serum (up to 1:1024) to determine the optimum concentration.
    - The host animal that produced the antiserum should also be run against the antiserum at its optimum concentration.
  - The antiserum is screened against a panel of common animal serums (diluted 1:100) to evaluate specificity. The host animal that produced the antiserum should also be run against the antiserum to evaluate non-specific cross-reactions.
    - If an antiserum reacts with non-related species or the host animal, the antiserum should not be used.
- 

**Quality control  
records**

The quality control tests for the normal serum and antiserum are recorded as follows:

- the dried gel is photographed and the image is imported into the *Crossover Electrophoresis Worksheet*.
  - the test is recorded along with the original gel, sealed in clear packaging, in the *Biology Quality Control Log Book* located in the Biology Laboratory.
- 

*Continued on next page*

## SER: Species Determination- Crossover Electrophoresis, Continued

### Sample preparation

Extract cutting, approximately 1 mm x 1 mm, of stained material in 50 µl deionized water, or place the cutting directly into the gel. Extracts containing blood should be a red-brown color.

### Procedure

Use the following procedure to perform crossover electrophoresis:

Step	Action
1	Pour buffer into tank. Allow buffer to reach room temperature.
2	Pour the support medium on the hydrophilic side of the gel bond, using a glass or plastic plate for support.
3	Punch wells.
4	Load 5 µl of antiserum into the anodic wells and 5 µl of sample, known standard, and blank into the cathodic wells.
5	Place the gel in the tank. Fill the tank with enough buffer to submerge the current lines on the gel and run the gel for a minimum of 30 minutes at 120 volts.
6	Remove gel from tank. Wash gel overnight in 1 M NaCl.
7	Rinse in deionized water, cover with filter paper (and a weight, if desired) to remove water, and dry in oven at approximately 60°C.
8	Place gel flat in a container and stain with <i>Stain Solution</i> for 2-5 minutes.
9	Destain the dried gel in <i>Destain Solution</i> until the background is clear and the precipitin bands are visible.
10	Photograph the stained gel and record the results on the <i>Crossover Electrophoresis Worksheet</i> . Import image into worksheet. Seal the original gel in clear packaging and store it with the evidence.

*Continued on next page*

## SER: Species Determination- Crossover Electrophoresis, Continued

---

**Interpretation** The following information must be considered when interpreting the results of this analysis.

- With a stereomicroscope, examine the area between the opposing wells. The presence of a stained precipitin band between two opposing wells is a positive (+) reaction.
    - Multiple precipitin bands may occur if the amounts of antigen and antibody are not balanced. Multiple bands do not negate a positive finding.
  - The absence of a precipitin band indicates a negative (-) reaction.
  - Cross-reactions with closely related species may occur. See product insert, located in the *Biology Quality Control Log Book*, for further information.
  - If bands are unclear or not in the proper position between wells, the test is inconclusive. The sample should be reanalyzed.
- 

**References** The following reference was used in the development of this procedure.

*Methods Manual*, Serological Research Institute.

---