

SER: Spermatozoa - Microscopic Examination

Principle

The cellular component of semen is spermatozoa. The microscopic identification of spermatozoa indicates the presence of semen.

Christmas Tree stain is a differential biological stain used to differentiate parts of cells and can be used to visualize spermatozoa. The two-reagent procedure stains nuclear material (sperm heads) red and epithelial membranes (sperm tails) green.

Equipment and supplies

This procedure uses the following laboratory equipment and supplies:

- compound microscope
 - centrifuge
 - oven/heat block
 - microscope slides
 - pipette
 - pipette tips
 - microcentrifuge tubes
 - ethanol
 - cover slips
 - mounting medium
-

Reagents

This procedure uses the following reagents:

- SERI Christmas Tree Stain Solution A
- SERI Christmas Tree Stain Solution B

Reagents can be stored at room temperature.

Record lot number and expiration date in examination records.

Continued on next page

SER: Spermatozoa - Microscopic Examination, Continued

Sample preparation

Use the following procedure to prepare the sample for slide examination.

Step	Action
1	Place one whole swab or stain cutting into a sterile microcentrifuge tube.
2	Add 300 µL of TE Buffer and incubate at room temperature for 30-60 minutes on the ThermoMixer.
3	Transfer 270 µl of the supernatant to a new tube.
4	Re-suspend the cellular pellet in the remaining 30 µl of supernatant.
5	Remove 2 µl of the re-suspended sample for slide examination.

Procedure

Use the following procedure to prepare and stain slides.

Step	Action
1	Dry 2 µl of re-suspended sample on a microscope slide for a minimum of 10 minutes in a 60° oven or on a heat block.
2	Stain with <i>Christmas Tree Stain A</i> for 5 to 30 minutes. Rinse with deionized water.
3	Stain with <i>Christmas Tree Stain B</i> for 5 to 60 seconds. Rinse with ethanol.
4	Allow the slide to dry. Add mounting medium and a coverslip.

Table continued on next page

SER: Spermatozoa - Microscopic Examination, Continued

Procedure
(continued)

Step	Action												
5	<p>Examine the slide with a microscope at 200x magnification. A phase contrast optical accessory can enhance visualization.</p> <ul style="list-style-type: none"> • With bright field microscopy, sperm heads appear red; the acrosomal region of the head appears lighter than other portions of the head. <ul style="list-style-type: none"> – Epithelial cells stain green, with the nuclei inside the epithelial cells appearing purple. • The sperm tail, if present, appears green. <ul style="list-style-type: none"> – With phase contrast microscopy, the acrosomal region appears darker than the rest of the sperm head. <p>NOTE: If the presence of debris or other cellular material interferes with the identification of spermatozoa, the analyst may clean up the sample and prepare a new slide.</p> <table> <tr> <th>Step</th><th>Action</th></tr> <tr> <td>1</td><td>Transfer 2 µl of the re-suspended sample to a new tube.</td></tr> <tr> <td>2</td><td>Add 300 µl of Digest Buffer and 15 µl of Proteinase K Solution.</td></tr> <tr> <td>3</td><td>Incubate for approximately 30 minutes at 56±1°C.</td></tr> <tr> <td>4</td><td>Pellet the cellular material by centrifugation.</td></tr> <tr> <td>5</td><td>Transfer the cellular pellet to a new slide, dry, and stain.</td></tr> </table>	Step	Action	1	Transfer 2 µl of the re-suspended sample to a new tube.	2	Add 300 µl of Digest Buffer and 15 µl of Proteinase K Solution.	3	Incubate for approximately 30 minutes at 56±1°C.	4	Pellet the cellular material by centrifugation.	5	Transfer the cellular pellet to a new slide, dry, and stain.
Step	Action												
1	Transfer 2 µl of the re-suspended sample to a new tube.												
2	Add 300 µl of Digest Buffer and 15 µl of Proteinase K Solution.												
3	Incubate for approximately 30 minutes at 56±1°C.												
4	Pellet the cellular material by centrifugation.												
5	Transfer the cellular pellet to a new slide, dry, and stain.												
6	Record how many sperm are present using the grading scale on the <i>Sexual Assault Examination Worksheet</i> .												