

## DNA: General Extraction Information

---

### Safety

### SAFETY WARNING

Biological samples have the potential to transmit infectious diseases. All analysts should use caution when performing this procedure and follow the *Bloodborne Pathogen Program* in the *Safety Manual*.

---

### Procedural safeguards

#### QAS 6.1.2

The following is a list of procedural safeguards for DNA extraction.

- Wipe down the work area with 10% bleach or ethanol before beginning the extraction process.
  - Put bench cover material (absorbent pads or butcher paper) on all work areas.
  - Never pipette reagents directly from a stock bottle. Reagents are to be poured into a sterile tube (for example, sterile disposable 50 mL centrifuge tube) before pipetting into sample tubes.
  - Have a wash bottle of 10% bleach solution or ethanol ready to wash all non-disposable implements (for example, forceps) each time they come in contact with samples. After washing with a stream of liquid, wipe dry with a clean lab-wipe tissue. Repeat twice.
  - Barrier pipette tips must be used when pipetting solutions containing DNA. Pipette tips must be changed between each sample, including the addition of the same reagent to each sample.
  - Wear gloves at all times and change them frequently when handling samples. Change gloves immediately if they come into contact with items containing DNA and may become contaminated.
  - Always extract reference samples separately from evidence samples.
  - Microcentrifuge tubes containing DNA should be pulse-centrifuged before removing the cap to minimize the possibility of creating DNA aerosols. Use a decapper to open tubes. Clean the decapper whenever it has come into contact with liquids containing DNA.
  - Do not touch the interior portion of the microcentrifuge cap.
- 

*Continued on next page*

## DNA: General Extraction Information, Continued

---

### Controls

The extraction procedures should use the following controls:

- QC stain control (with reference samples only)
  - reagent blank (at least two reagent blanks for all evidence sample batches)
  - substrate control (when appropriate)
- 

### Sample types

For extraction of each sample type, refer to the following documents:

- Non-hair samples such as blood, saliva, or contact samples
    - [DNA: Organic Extraction](#)
    - [DNA: BioRobot EZ1 and EZ1 Advanced XL Extraction](#)
  - Keratinized samples such as hair and fingernails
    - [DNA: Organic Extraction](#)
    - [DNA: BioRobot EZ1 and EZ1 Advanced XL Extraction](#)
  - Samples containing sperm
    - [DNA: Differential Extraction](#)
    - [DNA: BioRobot EZ1 and EZ1 Advanced XL Extraction](#)
    - [DNA: Organic Extraction](#)
-